

USERMANUAL

Chemical free UVC Room Disinfection

Disinfection step by step

1. Conduct the normal cleaning procedure
2. Make the room ready for disinfection
 - ✓ Move the beds, chairs and tables
 - ✓ Open the cabinet doors
 - ✓ Place the main unit in the center of the room
 - ✓ Place the satellites for irradiation of i.e., blind angles, toilets or entrance
 - ✓ Place the caution signs in front of the door/doors
3. Plug in the cable – earth connection
4. Push the green start/stop button
5. Choose time interval and push OK
6. Leave the room and close the door(s)
7. Check the display if the disinfection has been proceeded correctly – shown by ✓



NB!

You have 30 seconds to leave the room from you hit the start button until the disinfection process starts.

Disinfection times

Based on results from Department of Hygiene, Rigshospitalet, University Hospital of Copenhagen, we recommend:

- minutes exposure for non-spore-producing bacteria and viruses
- minutes for spore-producing bacteria and mycobacterium



SAFETY WARNINGS

The process can be stopped by

- Pushing "STOP UVC" at the caution sign
- Stopping the power outlet
- Pushing the green start/stop button at the main unit
- Pushing the red emergency stop on the top of the main unit
 - must be deactivated before next disinfection process can start



NB!

Use normal glasses/safety glasses if you need to enter the room while the disinfection is in progress

Maintenance

To keep the high and documented efficacy the
must be with similar
components.

Contact your local distributor of Dolphin Care
products to order similar components or sign up for a
service agreement.



Contact information

Head-office

Dolphin Care ApS
11, Unionsvej
DK-4600 Koege

Tel: +45 7174 1379

Mail: info@dolphincare.net

WWW: dolphincare.net

Local representative

TBN

Tel: +

Mail:

WWW:



We make high end disinfection products